

FOR BETTER CARE...

Jatab 'MM

THE MULTIFUNCTIONAL MOBILIZER



WITH THE PATIENT *and* THE CAREGIVERS IN FOCUS

Jatab 'MM is a unique assistive device for patient mobilization, transfer, and getting about.

The need for lifts and stressful tasks is minimized. Even toilet visits can be managed without moving the patient to another assistive device.

The risks of injuries from falls are reduced, as are the risks of occupational injuries.

Jatab 'MM can be adjusted to different heights, inclinations, and angles for all sitting positions from lying down to standing up. The chair is maneuvered electrically, and is very easy and safe to use.

JATAB CARE

When mobilization training is started early, the treatment results are strongly improved. The risks of complications are reduced, and the care period is shortened.

Jatab 'MM makes care tasks simpler and less arduous. Transfer can be done without lifting. The chair is continuously adjustable from lying to sitting. In case of blood pressure drop, it can be quickly adjusted to a lying position

With continuous position changes, the risks of pressure sores and hip contractures are reduced.

Mobilization



The importance of standing up



The standing-up function of the chair permits gradual training in a standing position.

This training is necessary for maintaining bowel activity, breathing, and blood circulation. It also counteracts spasticity.

In a standing position the patient can attain the same position as other standing people. This has great importance for emotional recuperation and self-confidence in interacting with others.

- Very easy to use
- Safe for the patient
- Easy on caregivers
- Dependable

- Better treatment results
- Shorter care period
- Lower cost

LESS LIFTING AND MORE CARE...

Jatab 'MM

MULTIFUNCTIONAL MOBILIZER

Transfer between bed and wheelchair is done with the patient in a prone position, using a special transfer board and sliding pad.

The transfer is accomplished with minimal strain and no lifting, and is very easy on the patient as well.

The working height is adjustable between 70 and 90 cm.

Transfer



Moving



Moving patients for operations, treatments, physiotherapy, dentist visits, etc. is accomplished smoothly in a **Jatab 'MM**, either lying down or sitting. With the sliding board the patient is easily transferred to a bed or an operation table.

Jatab 'MM has a built-in WC seat. Toilet visits can be simply and quickly accomplished. No lifts or transfers using other assistive devices are needed. This feature is especially important for visits outside the normal care environment.

Using the shower cover, the **Jatab 'MM** functions as a shower chair or shower gurney.

Hygiene



Jatab 'MM is a multi-functional assistive device that can handle all transfers, standing-up training, toilet visits, and showers without heavy lifts or lift devices.

Jatab 'MM features comprehensive treatment functionality that considerably reduces the risk of injuries caused by sitting too long in the same position.

Jatab 'MM offers comfort and a care environment permitting self-respect. The user can adjust the chair to his/her own comfort requirements

Jatab 'MM makes the caregiver's tasks easier, reduces heavy lifts, and thus reduces the risks of attritional injuries for the caregiver.

Jatab 'MM is available in a manual model or an electrically powered one.

Accessories

Upholstery	Removable, washable, optional incontinence protection
Sliding board	A special low-friction board
Sliding foot board	Facilitates transfer for a single caregiver
Sliding pad	Of washable fabric
Drain pan	Plastic drain pan with hanger mountings
Shower cover	With the shower cover, the chair can be used as a shower gurney
Adjustable table standing	Made of clear plastic with a spill protector. Can be used both sitting and standing
Electrical drive	Synchronized worm gear motor with joystick control

Powder-coated aluminum frame. Electric-powered actuators.
12V batteries with 220V battery charger. No-pump rubber wheels. Seat and back cushions made of polyurethane integral foam with nonporous surface. Removable, washable covers.
Propulsion motor 24V, maximum speed 6 km/hr. Maximum range about 20 km. (Extra device)

Jatab 'MM is CE-marked

Maximum width 680 mm
Seat width 500 mm
Seat depth 445 mm
Seat height 450 - 850 mm
Chair weight 48 kg
User weight 150 kg

JATAB CARE

Jatab Care AB • Drottninggatan 42 • Box 1540 SE-701 15 • Örebro • Sweden
Tel +46 (0)19 611 99 42 • Mob +46(0)70 557 25 58 • Fax +46(0)19 611 99 49 • Mail info@jatabcare.se

www.jatabcare.se